




Microbiome sample management

Collect, stabilize & extract DNA for gut microbiome profile studies

INVITEK
Molecular

Sample Collection Tubes with DNA Stabilizer



- ✓ Easily self-collect a high-quality sample using the pre-filled collection tubes
- ✓ Rapidly homogenize and stabilize at the point-of-collection
- ✓ Transport and store stabilized DNA at ambient temperature for 3 months – no cold chain required
- ✓ Eliminate the costs associated with temperature-controlled shipping
- ✓ Minimize bias introduced by microbial growth and DNA degradation

Product name	Package size	Catalogue number
Stool DNA Stabilizer	250 ml	1038111100
Stool Collection Tubes with DNA Stabilizer	50 tubes	1038111200
	250 tubes	1038111300

Extraction kits

Invitek Molecular is committed to developing products for state-of-the-art sample collection and nucleic acid extraction used in microbiome research. Maintaining the integrity of samples collected is a major challenge. Our collection and extraction systems guarantee an excellent performance in NGS downstream applications.




- ✓ Remove all PCR inhibitors from DNA samples
- ✓ Isolate DNA from microorganisms or from the host organism
- ✓ Extract nucleic acids from food and feed residues of plant or animal origin from the stool sample
- ✓ The kits are CE-marked and recommended for in-vitro diagnostic use
- ✓ Obtain high quality DNA suitable for 16S rRNA Microbiome profiling, Shotgun Metagenomic Sequencing, qPCR and arrays

Fast and convenient handling

The kits contain the InviAdsorb matrix in order to remove all effective inhibitors very efficiently. After a rigorous pre-lysis step using Zirconia Beads with optimized pre-lysis buffer at high temperatures, the sample is pre-incubated with InviAdsorb to remove PCR inhibitors. Proteinase K digestion ensures high yields also from Gram-positive bacteria. Stool samples contain a range of DNA e.g. host DNA from colon epithelial cells, parasite DNA, bacterial DNA, DNA from food or DNA from gastrointestinal pathogens. The choice of different lysis conditions enables the accumulation or reduction of the bacterial DNA content in the total DNA.

DNA extraction after sample storage at ambient temperatures



1.4 ml of preserved stool samples were tested on day 1, day 2, day 7, day 18, day 49 and day 105 for DNA isolation using the PSP Spin Stool DNA Plus Kit. 2 µl eluate were used in the RIDA®GENE Gut Balance real-time PCR assay from R-Biopharm

Kit	Starting material	Package sizes	Catalogue number
Membrane adsorption			
PSP® Spin Stool DNA Kit	up to 200 mg	50 preps	1038100200
		250 preps	1038100300
PSP® Spin Stool DNA Plus Kit	1.4 ml Stool DNA Stabilizer with stool homogenate	50 preps	1038110200
		250 preps	1038110300
PSP® Spin Stool DNA Basic Kit (without Stool Collection Tubes)	1.4 ml Stool DNA Stabilizer with stool homogenate	50 preps	1038120200
		250 preps	1038120300
Automated using magnetic beads			
InviMag® Stool DNA Kit/KF96 (for use on KingFisher™ 96 and KingFisher™ Flex, Thermo Fisher Scientific)	up to 200 - 400 mg	5x96 preps	7438300200



Compliance with EU Directive 98/79/EC on in vitro medical devices. Not for in-vitro diagnostic use in countries where the EU Directive 98/79/EC on in vitro medical devices is not recognized.

PSP® and InviMag® are registered trademarks of Invitek Molecular GmbH. Other brands or product names are trademarks of their respective holders. © 2020 Invitek Molecular, all rights reserved.

6If/09/2020

INVITEK
Molecular

Invitek Molecular GmbH
Robert-Rössle-Str. 10
13125 Berlin Germany

Phone: +49 30 9489 2908
info@invitek-molecular.com
www.invitek-molecular.com